



DoverPac® System Compared to a "Big Bag"



A proven, complete system is critical in your HAZOP review to eliminate risk of failure.

DoverPac®	vs.	Common Big Bag
▶ Contained Attachment for Filling & Emptying		▶ Not Designed for Contained Applications
▶ Contained Tamper Proof Sampling		▶ Must Open for Sampling
▶ Low Capital Cost		▶ Requires Expensive, Maintenance Intensive Capital for Contained Attachment
▶ 100% Inflation Tested and Inspected		▶ No Inflation Acceptance Testing and Limited Inspection...Not High Containment Safe
▶ Custom or Standard Sizes		▶ Standard Sizes
▶ Contained Separation Using Crimps		▶ Wire Ties Only or Twist, Tie, Tape, Cut Separation
▶ System Application Results in Validated Containment Levels		▶ No Long Term Verification in Contained Applications
▶ Type C Restraint made with rugged ArmorFlex® Film		▶ Type C Restraint
▶ CE Mark on liner and restraint		▶ No CE Mark on liner
▶ Audited Quality and Manufacturing Systems		

Beyond Boundaries™.

Innovators at our core, we develop engineered solutions for our customers' complex problems. Recognized globally for our flexible containment solutions, ILC Dover serves customers in a diverse range of industries, including pharmaceutical and biopharmaceutical manufacturing, personal care, food and beverage, chemical, aerospace, healthcare and government agencies. At ILC Dover, quality is a culture, not a measurement. Our customers will tell you that we cater to their every need and that we're highly innovative, responsive, dedicated and competitive. We have been innovating since 1947. ILC Dover's visionary solutions improve efficiency, safeguard workers and product, and prevent disasters – proof that we are on the front line of business excellence.

Engineering evolution beyond boundaries.



One Moonwalker Road
Frederica, DE 19946 USA
+1.302.335.3911
+1.800.631.9567

customer_service@ilcdover.com

www.ilcdover.com

