

Flexible Containment Solutions Guide

ILC DOVER
creating what's next»

DoverPac
containment solutions»


Eclipsing the traditional FIBC, the G2Pac® is engineered to address the issues that arise from uncontrolled processes, including compromised facility cleanliness, cross contamination risks and employee exposure to unnecessary hazards. G2Pac® supports cGMP operations.

» FEATURES

- Durable ArmorFlex® liner with groundable restraint (type C)
- Individual necks with bag-out sleeve to support charging and offloading operations
- Secure, dust-free hardware interface to process equipment
- Sized to interface with standard pallet
- Interfaces with ILC contained separation crimps

» BENEFITS

- Economical containment
- Supports dust-free operation
- Reduced cross contamination risks
- Protects operator
- One-person operation



| Feature | |
|---|------------------|
| Restraint and Liner | Yes |
| Charge/Offloading Neck | Yes |
| ArmorFlex® 104 | Yes |
| Sizing | 700L, 560L, 185L |
| Integral Bag-Out Sleeves for Contained Changeover | Yes |
| Interfaces with Contained Separation Crimps | Yes |

MKT0065 Rev C

