



Sentinel XL[®] HP



HP SYSTEM XL BLOWER — ROBUST PROTECTION AGAINST A RANGE OF THREATS



BENEFITS YOU CAN EXPECT

Designed to protect healthcare workers from diseases/
pandemic flu.

- Single size head cover eliminate need to buy and stock two sizes.
- Head cover design allows stethoscope use.
- Audible/visual low-voltage alarm.
- NiMH batteries eliminate “memory effect” associated with NiCd batteries.
- Rugged dual chamber blower design for increased durability.
- Secure threaded hose/headcover interface is easy to use.
- NIOSH approved



One Moonwalker Road
Frederica, DE 19946 USA
+1.302.335.3911
+1.800.631.9567
customer_service@ilcdover.com

www.ilcdover.com