



# Flexible Containment Solutions Guide



Eclipsing the traditional FIBC, the G2Pac® is engineered to address the issues that arise from uncontrolled processes, including compromised facility cleanliness, cross contamination risks and employee exposure to unnecessary hazards. G2Pac® supports cGMP operations.

## FEATURES

- Durable ArmorFlex® liner with groundable restraint (type C)
- Individual necks with bag-out sleeve to support charging and offloading operations
- Secure, dust-free hardware interface to process equipment
- Sized to interface with standard pallet
- Interfaces with ILC contained separation crimps

## BENEFITS

- Economical containment
- Supports dust-free operation
- Reduced cross contamination risks
- Protects operator
- One-person operation

Feature	
Restraint and Liner	Yes
Charge/Offloading Neck	Yes
ArmorFlex® 104	Yes
Sizing	700L, 560L, 185L
Integral Bag-Out Sleeves for Contained Changeover	Yes
Interfaces with Contained Separation Crimps	Yes

One Moonwalker Road  
Frederica, DE 19946 USA  
+1.302.335.3911  
+1.800.631.9567  
customer\_service@ilcdover.com

[www.ilcdover.com](http://www.ilcdover.com)