

Sentinel EZ BioHood™

Protective PAPR Hood for Healthcare Workers on the Frontlines

COVID-19 and other infectious diseases threaten the safety of heroic healthcare workers. The Sentinel EZ BioHood by ILC Dover offers 100x the protection of an N-95 mask, providing superior respiratory protection to healthcare workers on the frontline and enabling them to provide continued care to patients. The hood's unique design makes it intuitive to use, easy to use, store, and workers do not need to be fit-tested to wear it effectively.

FEATURES

Superior protection

- 1,000 APF* provides 100x the protection of an N-95

Readily Available

- Produced in our Frederica, DE facility, this hood is available for rapid shipment.

No fit-testing required for wearers

- Adjustable up to 5 sizes in 1 single-size product
- Reduces logistic burden of carrying multiple sizes

Easy to store

- Compact, lay flat design frees storage space by 50%

Comfortable to wear for long operations

- Breathable loose-fitting neck cuff
- Internal head band provides good stability
- Provides ability to use stethoscope

Wide range of view

- Features a wraparound Clear Visor for superior visibility
- Healthcare worker is clearly visible to the patient

Infectious Disease PAPR

Sentinel EZ BioHood is NIOSH approved (TC-21C-1130) as a loose-fitting hood designed to work with the Sentinel XL® blower system.



Flexible Protective Solutions

We improve efficiency while safeguarding people, product, and infrastructure in hazardous conditions.

LET'S GO
BEYOND BOUNDARIES
TOGETHER.

One Moonwalker Road
Frederica, DE 19946 USA
+1.302.335.3911
+1.800.631.9567
customer_service@ilcdover.com

www.ilcdover.com

*Based on manufacturer data