



## Grayling™ Liners for Steel, Plastic and Fiber Drums

### **GRAYLING™ DRUM LINERS FROM ILC DOVER**

An Efficient, Economical Solution to the Cleaning/Disposal Process  
for Used Drums

Grayling drum liners are manufactured using engineered polyethylene film structures to provide durability, consistency and performance and are manufactured in an ultra clean environment to ensure that your product stays pure and contamination free.

ILC Dover keeps a large inventory of standard sized drum liners in stock and we will also custom manufacture liners for your unique application.

ILC Dover is a leading manufacturer of IBC liners for industrial processors of bulk materials. You can depend on ILC Dover's Grayling liners to protect your product during processing, storage and equipment.

# Grayling™ Liners for Steel, Plastic and Fiber Drums

## STOCK LINER SPECIFICATIONS

PART #	DESCRIPTION	DIMENSIONS (IN)	GAUGE	CONSTRUCTION/ MATERIAL TYPE	QTY
920344	55 gallon round bottom cuff top drum liner	23" DIA x 40" tall	4MIL	LDPE TYPE 43	75/CASE
920345	55 gallon round bottom twist top drum liner	23" DIA X 56" TALL	4MIL	LDPE TYPE 43	50/CASE
920346	55 gallon round bottom cuff top drum liner	23" DIA X 40" TALL	8MIL	LDPE TYPE 43	40/CASE
920347	55 gallon round bottom twist top drum liner	23" DIA X 56" TALL	8MIL	LDPE TYPE 43	30/CASE
920348	55 gallon round bottom twist top drum liner	23" DIA X 56" TALL	3MIL	HIGH HEAT PP TYPE 30	75/CASE
921660	55 gallon round bottom twist top drum liner	23" DIA X 56" TALL	3.5MIL	BARRIER FOIL TYPE 88	50/CASE

SPECIAL APPLICATION LINERS ARE ALSO AVAILABLE:

- Anti-Static Liners
- Chemical Resistant Liners (Nylon)

Beyond Boundaries.

Inventors at our core, we develop engineered solutions for our customers' complex problems. Recognized globally for our flexible containment solutions, ILC Dover serves a diverse range of customers, including pharmaceutical and biopharmaceutical manufacturing, personal care, food and beverage, chemical, aerospace, healthcare, and government agencies. At ILC Dover, quality is a culture, not a measurement. Our customers will tell you we cater to their every need and that we are highly innovative, responsive, dedicated and competitive. We have been inventors since 1947.

ILC Dover visionary solutions improve efficiency, safeguard workers and product, and prevent disasters — proof that we are on the frontline of business excellence. Engineering evolution beyond boundaries.



ILC Dover LP  
100 Discovery Blvd, Suite 1001  
Newark, DE 19713 USA  
+1.302.335.3911  
+1.800.631.9567  
customer\_service@ilcdover.com

[www.ilcdover.com](http://www.ilcdover.com)