soloFLEX Benchtop Containment Isolator



High Performance. Rapid Assembly. Fast Delivery.

Meet the soloFLEX benchtop containment isolator. Complete with 3x H14 HEPA filters, material entry airlock, and a continuous liner port, this standard isolator system is capable of handling most lab scale processes. It is available in three standard formats: 2 Glove, 3 Glove, and 4 Glove, and can be easily assembled in under 30 minutes. Order the full isolator system or a replacement enclosure today.



Small & Compact Benchtop Isolator

The complete isolator includes flex enclosure and compact exhaust fan. The enclosure is disposed of as a complete unit using the fan vacuum to collapse. Replacement enclosure systems are supplied as a ready to assemble kit.

- Robust stainless steel support frame
- 3 x H14 HEPA filters
- 3x fan options
- 300mm (12") Dia. continuous liner x 130Ft liner cassette
- Manual and automatic speed control fans with 120V/60 supply

- Compressed air venturi fan available
- Low friction polyurethane glove sleeves
- Utility sleeves for utility cables 2" Dia.
- PVC base plates for balance stability

SAFE, EASY, CONVENIENT

Rapid and easy assembly

Modular stainless frame

High containment performance

Contamination-free disposal

Low price and fast delivery



How soloFLEX Works

High-Tension Flexible Film Enclosure



soloFLEX operates with a high-tension flexible film enclosure that permits operation at significant negative pressures, preventing hazardous API materials from escaping the contained space.

Double-Zipper Airlock Protection



soloFLEX features double-zipper airlock protection for entering materials into an ongoing process. Using proven SOP's, soloFLEX is easy and intuitive and offers excellent containment performance.

Convenient & Captive HEPA Filters



soloFLEX's technical components are bonded to the flexible enclosure, allowing for an easy, contamination-free disposal of the contaminated enclosure, without the need for any cleaning.

Matereials / Operations	Yes	No
Pharma powder sampling	✓	
Pharma weighing and dispensing	~	
Pharma IPC Instrument housing	✓	
Small process equipment housing	✓	
Cytotoxic materials		Factory Special
Organic solvents		×
Hot processes/materials		×
Inert gas		Factory Special
Drums up to 12" Dia. x 12" H	~	
Drums over 12"		×



Low Price, Fast Delivery! Contact Solo Containment for Ordering Information

Buy the complete isolator system or a replacement enclosure system.



Manual Speed Control 120V/60Hz 1Ph x 70w



Automatic Speed Control 120V/60Hz 1Ph x 70w



Manual Needle Valve controlled compressed air venturi





Solo Containment is now an **ILC Dover LP Company**

Technology Park Unit 3 Bramhall Moor Ln Pepper Rd, Hazel Grove Stockport SK7 5BW United Kingdom +44(0)161 358 2233

www.solocontainment.com

sales@solocontainment.com