



# ArmorFlex® 114 Film

At ILC Dover we are going “beyond boundaries” by listening carefully to our customers. In response to customer input and emerging regulations, we’ve developed a robust new film to support our DoverPac® Containment Systems. This revolutionary film, ArmorFlex® 114, delivers permanent static dissipative properties while meeting FDA and 2002/72/EC requirements for food contact and EP 3.1.3 test conditions.

- Complies with FDA 21 CFR
- Fully compliant to EU Regulation 10/2011 and amendments
- Meets test parameters of EP 3.1.3
- Passes USP <661> Physicochemical Tests for Plastics
- Passes USP <88> Class VI (7 day implant)
- Passes <87> biological reactivity, in vitro
- 5-year shelf life
- Permanent antistat replaces migrating additives
- Passes Chilworth incendivity tests
- Tested for solvent resistance
- DMF filed with FDA



SCAN ME



BEYOND BOUNDARIES™

ILC Dover LP  
One Moonwalker Road  
Frederica, DE 19946 USA  
+1.302.335.3911  
+1.800.631.9567  
customer\_service@ilcdover.com

[www.ilcdover.com](http://www.ilcdover.com)



MKT-0049 Rev C

\*ArmorFlex and DoverPac are registered trademarks of ILC Dover LP.